


biocomma® PCR Plate Sealing Films



biocomma® PCR Plate Sealing Films are used for sealing high throughput plates, which also optically clear for qPCR tests. Adhesive films and pressure-sensitive films are available to meet different needs.

Features:

- Compatible with 96 and 384 well plates
- Superior good sealing of each well
- High fluorescence transmittance across common fluorophores
- Minimize evaporative loss
- Strong adhesive, suitable for heat sealer
- Easy to apply manually
- DNase & RNase Free



Ordering Information

Cat. #	Description	Qty.
96WT-DC	biocomma® PCR Plate Sealing Films, strong adhesive, highly optically clear, two easy removal white end tabs, 137*80 mm	100 Pcs/PK
96WT-DC-1	biocomma® PCR Plate Sealing Films, pressure-sensitive, highly optically clear, two easy removal white end tabs, 140*80mm	100 Pcs/PK

Please contact us for more specifications.