



biocomma



## Environmental Sampling Kit

# Up to 95% sample collection and release

### Features:

It is used for sampling on the surface of equipment in public health places, equipment in hospitals, equipment in cosmetics industry, equipment in food industry, and its operating table surface and operators' hands. And it also can be used for sampling on the surfaces of public health areas where is containing preservatives and disinfectants.

The main ingredients of CY001E are Sodium Chloride and distilled water. It is used for sampling on the surface of unsterilized objects.

The main ingredients of CY002E are Sodium Chloride, Sodium Thiosulfate, Tween 80, Lecithin, and glycine compound neutralizer. It is used for sampling on the surface of sterile objects (within 2 hours).



### Order Information:

Cat.#	Product Name	Description	Qty.
CY001E	Environmental Sampling Kit	10mL preservation tube with 3mL preservation medium, 1* sampling swab	50 pcs/box
CY002E	Environmental Sampling Kit (including neutralizer)	10mL preservation tube with 3mL preservation medium (containing neutralizer), 1* sampling swab	50 pcs/box

# BIOCOMMA LIMITED

Add: Ground FL, Bdg 12,Zhonghaixin Industrial Park ,Ganli Six Rd, Jihua St, Long Gang Dist , ShenZhen, 518114 P.R.China.

TEL:86(755)-25431879

FAX:86(755)-25498726

WEB:www.biocomma.com

Email:info@biocomma.com