


Virus DNA/RNA Extraction Kits (Pre-Packed Magnetic Beads)

Great Upgrade!



**No Need of
Proteinase K Anymore**



Efficiency

Comparing to traditional extraction method, the time of dispensing and adding Proteinase K to each well are saved. The whole extraction process could accelerate **30~45%**.

Economy

Since no refrigeration needed for storage and transportation anymore, the cost is significantly reduced, and storage space is also saved a lot.

Reproducibility

As shown on figure of Real-time PCR reaction, the curves without Proteinase K express the same as traditional method. All results of Ct value are strong positive, so that let you test with confidence.