

Manual Pipettes

- Wide measuring range from 0.1µl to 10ml.
- EN/IOS08655 standard calibration, and online calibration is available.
- The lower half can be sterilized at high temperature and pressure.

Ordering Information:

Cat. #	Channels	Range
7010101001	1	0.1-2.5µl
7010101004	1	0.5-10µl
7010101006	1	5-50µl
7010101005	1	2-20µl
7010101008	1	10-100µl
7010101009	1	20-200µl
7010101014	1	100-1000µl
7010101016	1	200-1000µl
7010101017	1	1000-5000µl
7010101033	1	2-10ml
7010103004	8	0.5-10µl
7010103006	8	5-50µl
7010103012	8	50-300µl



Bottletop Dispensers

- Chemical resistance, meeting the capacity range of 0.5-100mL.
- Various accessories are available; partial or whole sterilization at high temperature is available.
- Backflow valve design.

Ordering Information:

Cat. #	Range
7032100001	0.5-5ml
7032100002	1.0-10.0ml
7032100003	2.5-25.0ml
7032100004	5.0-50.0ml
7032212005	10-100ml
7032212004	5.0-50ml
7032212003	2.5-25ml
7032212002	1.0-10ml
7032212001	0.5-5ml

