




CommaClean[®] Square Media Bottles

Features:

- Made of, light weight, and anti-crushing PET
- Highly transparent and thickened bottle body is easy to observe
- Good chemical stability is suitable for storing high-clean reagents, granular standard substances and standards
- Founded bottle body conforms to ergonomics, and easy to hold
- Unique T-shaped thread structure to ensure good sealing
- Strict leak test to ensure no leakage during air transportation
- Produced in Class 100,000 GMP workshop
- Gamma ray sterilization, up to 10⁻⁶ Sterility Assurance Level (SAL)
- Strict quality control on raw material and products to meet biological testing standards
- Used to pack and transport liquid medium, buffer and serum
- Pyrogen free





Applications

CommaClean® Square Media Bottles are suitable for storage or transportation of reagents, especially for production and storage of cell culture medium, cell serum, sterile buffer and other reagents.



All caps are universal. Also, universal for the most manufacturers in the market.



Injection, stretch and blow molding at one time to ensure good sealing.



Ergonomic founded design is easy to hold and save space.



Good consistency of products.

Cat. #	Volume (mL)	Length (mm)	Width (mm)	Height (mm)	Diameter (mm)	Weight (g)	Qty.
22103	50	40	40	82	22	18	72 Pcs/PK 5 PK/Ctn.
22203	125	54	54	108	28	46	25 Pcs/PK 4 PK/Ctn.
22303	250	59	59	146	28	54	30 Pcs/PK 2 PK/Ctn.
22401	500	76	76	175	28	87	24 Pcs/PK 2 PK/Ctn.
22503	1000	92	92	216.5	28	114	24 Pcs/PK 2 PK/Ctn.