


DC10H

Automated Screw Cap Capper/Decapper



Compatibility

Compatible with multiple specifications of collection tubes

Intelligent

The opening speed is adjustable
Tube position is detected by infrared sensor
Automatically record the number of operations

Efficient

One-hand operation
Complete capping or decapping within 3s

Portable

Dimension is only 18*22*50 cm
2 units can be placed in 1 biological safety cabinet

Cat. #	Description	Qty.
BM-DC10H-E	Automated Screw Cap Capper/Decapper	1 pc/box