

Features

- USP Class VI and FDA certified medical grade polypropylene raw materials.
- Thin tube wall, fast heat transfer, precise temperature control.
- Ultra-transparent tube cover design can ensure that the fluorescent signal can pass through without obstacles.
- Good sealing performance of the tube cover, which prevents contamination, and is easy to open.
- Autoclavable (121°C sterilized for 20min).
- DNA-free, DNase-free, RNase-free and pyrogen-free.

PCR Tube/Plate

Efficient

Precise

Ordering Information

Cat.#	Description	Qty.
ST002	0.2 mL, single PCR tube, transparent, flat cap	1000pcs/bag
ST005	0.5 mL, single PCR tube, transparent, flat cap	1000pcs/bag
4WT-001	0.1 mL, 4-strip PCR tube, transparent	250sets/box
8WT-001-CHN	0.1 mL, 8-strip PCR tube, transparent, suitable for 8-strip flat cap	125sets/box
8WT-002-CHN	0.2 mL, 8-strip PCR tube, transparent, suitable for 8-strip flat cap	125sets/box
8WT-001-1	0.1 mL, 8-strip PCR tube, transparent, suitable for 8-strip dome cap	125sets/box
8WT-002-1	0.2 mL, 8-strip PCR tube, transparent, suitable for 8-strip dome cap	125sets/box
8WT-001-2	0.1 mL, 8-strip PCR tube with attached individual caps, transparent	125sets/box
8WT-002-2	0.2 mL, 8-strip PCR tube with attached individual caps, transparent	125sets/box
8WT-001W	0.1 mL, 8-strip PCR tube, white, suitable for 8-strip flat cap	125sets/box
8WT-001W-1	0.1 mL, 8-strip PCR tube, white, suitable for 8-strip dome cap	125sets/box
8WT-002W-1	0.2 mL, 8-strip PCR tube, white, suitable for 8-strip flat cap	125sets/box
8WT-002W-2	0.2 mL, 8-strip PCR tube, white, suitable for 8-strip dome cap	125sets/box
96WT-001-N	0.1 mL, 96 well PCR Plate, transparent, non-skirt	25pcs/box
96WT-001-H	0.1 mL, 96 Well PCR Plate, transparent, semi-skirt	25pcs/box
96WT-001W-N	0.1 mL, 96 Well PCR Plate, white, non-skirt	25pcs/box
96WT-001W-H	0.1 mL, 96 Well PCR Plate, white, semi-skirt	25pcs/box
96WT-002-N	0.2 mL, 96 Well PCR Plate, transparent, non-skirt	25pcs/box
96WT-002-H	0.2 mL, 96 Well PCR Plate, transparent, semi-skirt	25pcs/box
96WT-002W-N	0.2 mL, 96 Well PCR Plate, white, non-skirt	25pcs/box
96WT-002W-H	0.2 mL, 96 Well PCR Plate, white, semi-skirt	25pcs/box
384WT-001-H	30 μ L, 384 Well PCR Plate, transparent, full-skirt	25pcs/box
384WT-001W-H	30 μ L, 384 Well PCR Plate, white, full-skirt	25pcs/box

